



CLEAR ALIGNER

CA® TIP Pliers

CA® TIP 1

Large spherical tip allows insertion of retention points in CA® splints, for increase of friction. Stainless steel, sterilisable.
Warranty: 5 years.

CA® TIP 1	5281.1 1 piece
-----------	-------------------

CE

CA® TIP 2

Small spherical tip allows forming of activation points in CA® splints. Stainless steel, sterilisable.
Warranty: 5 years.

CA® TIP 2	5291.1 1 piece
-----------	-------------------

CE

CA® TIP 3

Pliers with oval plateau, enable creation of activation surfaces in CA® splints. Stainless steel, sterilisable.
Warranty: 5 years.

CA® TIP 3	5292.1 1 piece
-----------	-------------------

CE

CA® TIP 4

Pliers for torque control with curved plateau, enable creation of labial/lingual activation surfaces in CA® splints. Stainless steel, sterilisable. Warranty: 5 years.

CA® TIP 4	5295.1 1 piece
-----------	-------------------

CE

CA® TIP 5

Pliers for aligners with angled pliers head for correction of splint edges in case of gingival pressure marks. Stainless steel, sterilizable.
Warranty: 5 years.

CA® TIP 5	5299.1 1 piece
-----------	-------------------

CE

CA® TIP 6

Pliers for aligners for removal of activation points and surfaces inserted by means of CA® TIP 2, 3 and 4. Also suitable for finishing and smoothing of the splint surface. Stainless steel, sterilizable.
Warranty: 5 years.

CA® TIP 6	5323.1 1 piece
-----------	-------------------

CE

